

Controlled Drugs Kit (CD Kit)



The PHS Wastemanagement Controlled Drug (CD) Destruction Kit is recognised as the securest way to denature drugs to render them safe and without value.

By placing such waste in these containers, you can rest assured that it is stored correctly and safely.

Capacities & Dimensions

250ml	11.2cm x 6.4cm dia
2l	20.7cm x 11.8cm dia
4l	24.8cm x 15.5cm dia
20l	43.8cm x 24.5cm dia

Instructions for use:

Controlled drugs should be removed from blister packs where necessary and crushed. Crushed tablets and any liquids should be emptied into the disposal kit, and water added as indicated.

The contents should then be shaken until fully mixed. The whole kit should subsequently be disposed of in the relevant waste bin.

Health and Safety Precautions: Use these containers as per the instructions, keeping all waste securely within the container.